) **Velat[®] Ge**

mucopolysaccharide, adrenocortical ext., salicylic acid

cooling, soothing & effective



- Direct to the pain
- Direct to the swelling

Movelat Gel

The new 100g tube adds both economy and convenience.

- Tried
- Tested
- Trusted

The Direct Anti-inflammatory

Presentation: Movelet Gel contains mucopolysecheride polysulphuric acid ester (heperinoid) 0,2g, adrenococtical extract 1,0g (equivelent to 20mg of corticosteroids), and salicylic acid 2,0g in a colourless gel base to 100g. Movelet Cream is also available containing the earne active constituents in a

Indications: Non-rheumatic disorders of the locomotor apparatus: sprains, muscular strains and spasms, soft tissue injuries.
Rheumatic disorders: osteo-arthritis, fibrositis.

Dosege and administration: Movelat Gel: Two to six inches (5-15cm) to be applied gently to the affected area up to four times daily; may also be used as an ultrasound coupling medium. Movelat Cream: Two to six inches (5-15cm) to be massaged into the affecte area up to four times daily; may also be used as a dressing, when it should be applied thickly.

Contra-Indications, warnings: Movelet Gel: Contains an alcohol and therefore should not be brought into direct contact with open wounds or mucous membranes Movelat Cream: Nil

£1,97 per 50g tube of Gel or Cream £3,90 per 100g tube of Gel or Cream

0542/0006 Holders and Distribu

Full information on request: LUITPOLD-WERK Medical & Scientific Office in U.K. Hayes Gete House, 27 Uxbridge Road, Hayes, Middlesex UB4 0JN, Tel: 01-561 8774

